



Disinfecting the diagnostic monitors

The following suggested procedures can keep the products clean and hygienic.

Applicable products

Diagnostic and clinical monitors

Items to prepare

1. Mask, disposable gloves and protective glasses
2. Cleaning cloth with soft and lint-free materials such as microfiber cloth
3. Disinfectant/cleaning solution

Use one of the following is considered to effectively clean and keep the surface hygienic.

- a. 70% isopropyl alcohol
- b. 80% alcohol (ethyl alcohol)

Note 1: Do not use thinners, solvents, waxes or polishing cleaners.

Note 2: Refer to the manufacturer's instructions for safety use of the disinfectant/cleaning solution.

Disinfection procedure

1. Wear a mask, disposable gloves and protective glasses.
2. Turn off the product and disconnect the power cord.
Working on the product with the power supply connected may result in an electric shock or malfunctions.
3. Remove any other connected cables and equipment.
4. Soak the cleaning cloth with the disinfectant/cleaning solution.
Note) Do not let the liquid drip. Do not spray directly on the product.
5. Wipe the product gently and thoroughly. When working on the screen, do not push or rub it to avoid screen failures and/or scratches.
6. When completed, let the product dry completely.
7. Follow the above procedures to wipe the power cord, AC adapter and connection cables. If you need to work on other connected devices, follow each manufacturer's instructions.
8. Once finished, discard the mask, disposable gloves and protective glasses, then wash your hands.

Before turning on the product

1. Make sure the product, power cord and cables are completely dry.
2. Connect all the cables and equipment to the product.
3. Connect the power cord.

Improper work or misuse of a cleaning agent may lead to changes in product gloss, discoloration, fading, and/or deterioration of image quality. Experiment on an inconspicuous part beforehand would be highly recommended.

※Please [contact us](#) if you have any questions.